

# LLG-Carboy, PP, with handles

LLG-Carboy with screw cap made of PP for storing and dispensing media, solutions and other liquids. The practical shoulder handle facilitates transport and pouring.

- With graduation
- Autoclavable
- Stable bottom with curved corners
- Leak-proof

Cap size: 83 mm  
No. of handles: 2  
With stopcock: No

Nominal capacity l	Ø mm	Height mm	Dimensions (ØxH)	Neck int. Ø mm	Colour	Material	Autoclavable	PK	Cat. No.
10	249	376	249 x 376 mm	65	Natural	PP	Yes	1	4.692 534
20	285	518	285 x 518 mm	65	Natural	PP	Yes	1	4.692 535
50	376	665	376 x 665 mm	65	Natural	PP	Yes	1	4.692 536