

# Qualitative filter paper, Grade 593 $\frac{1}{2}$ ;, folded filters

## Medium to slow 2 to 5 $\mu\text{m}$

- time saving, because folded filters have a larger surface area than comparable circular filters
- for medium to fine precipitates and general laboratory work
- Weight: 85g/m<sup>2</sup>;
- Thickness: 0.17mm

Type	Ø mm	Weight g/m <sup>2</sup>	Thickness mm	Particle retention $\mu\text{m}$	Type of analysis	Shape	PK	Cat. No.
593 $\frac{1}{2}$ ;	185.0	85	0.17	2...5	Qualitative	Folded filters	100	9.063 106
593 $\frac{1}{2}$ ;	240.0	85	0.17	2...5	Qualitative	Folded filters	100	9.063 107

## Alternativen:

Filter papers 593, qualitative, folded filters / Hahnemühle FineArt