

Graduated jugs, SAN, moulded graduations

Crystal clear. With a raised scale and stable, easy-grip handle. Food-safe product according to regulation (EU) No. 10/2011.

Colour:

Crystal clear

Nominal capacity ml	Graduation ml	Ø mm	Height mm	With spout	With handle	Material	Type of graduation	Food safe	PK	Cat. No.
250	5	70	120	Yes	Yes	SAN	Embossed	Yes	1	7.008 966
500	10	91	133	Yes	Yes	SAN	Embossed	Yes	1	7.008 968
1000	10	116	170	Yes	Yes	SAN	Embossed	Yes	1	7.008 970
2000	20	150	215	Yes	Yes	SAN	Embossed	Yes	1	7.008 972
3000	50	170	242	Yes	Yes	SAN	Embossed	Yes	1	7.008 974

Alternativen:

Measuring jugs with handle, PP / Kartell S.p.A.