

Graduated jugs, PP, moulded, blue graduations

Highly transparent. With moulded, blue graduations, spout and stable, easy-grip handle.
 Food-safe product according to regulation (EU) No. 10/2011. To preserve markings, do not clean at temperatures exceeding 60 °C. Conditionally autoclavable at 121 °C (2 bar) according to DIN EN 285. For autoclaving, we recommend the design with raised graduations.

Colour:

Transparent

Nominal capacity ml	Graduation ml	Ø mm	Height mm	Standard	With rim	With spout	With handle	Material	Type of graduation	Color
50	2	40	70	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue
100	2	50	80	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue
250	5	74	120	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue
500	10	92	140	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue
1000	10	117	181	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue
2000	20	152	213	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue
3000	50	172	242	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue
5000	50	204	270	DIN EN 285	Yes	Yes	Yes	PP	Embossed	Blue

Alternativen:

Graduated jugs, PP, coloured / VITLAB

Measuring jugs with handle, PP / Kartell S.p.A.

Graduated jugs, PP, moulded graduations / VITLAB

Graduated jugs, PP, moulded graduations / VITLAB