

Ultrasonic bath BactoSonic^{®reg} 14.2 Set

Ultrasonic special bath for a gentle removal of biofilms.

Scope of supply: Ultrasonic bath, scientifically proven description of the procedure, transport containers PP: 2 x IB 5 (0.52 l), 2 x IB 6 (0.6 l), 1 x IB 10 (1.0 l), 1 x IB 18 (1.8 l), 1 x IB 20 (2.0 l), container carriers for transport containers PC: 2 x BT 5 for IB 5, 2 x BT 6 for IB 6, 1 x BT 10 for IB 10, 1 x BT 18 for IB 18, utensil holder for transport containers GH 14 (stainless steel, for 3 x IB 20), frame for foil test FT 14, 250 ml contact liquid TICKOPUR R 33 (concentrate)

Int. Dim. (WxDxH):	325 x 300 x 150 mm
Ext. dim. (WxDxH):	355 x 325 x 305 mm
With stopcock:	Yes
Weight:	10.2 kg
Power supply:	230 V, 50/60 Hz
IP code:	IP 32

Type	Nominal capacity l	Ultrasonic frequency kHz	Ultrasonic power W	Display	Type of display	Open
BactoSonic ^{®reg} 14.2	13.50	40	200	Analogue	Without display	Rotat
BactoSonic ^{®reg} 14.2	13.50	40	200	Analogue	Without display	Rotat

Zubehör:

Transport containers, PP for BactoSonic[®] 14.2 / SONOCOOL[®] 255 / Bandelin electronic

Container holders, PC for BactoSonic[®] 14.2 / Bandelin electronic

Concentrates for ultrasonic baths TICKOPUR R 33 / Bandelin electronic

Alternativen:

Ultrasonic bath BactoSonic[®] 14.2 / Bandelin electronic