

Bunsen burners with air regulator and needle valve

With air regulator and needle gas valve. Adjustable nozzle.

Combustible	Height mm	With air regulation	Weight g	PK	Cat. No.
All-gas	170	Yes	270	1	9.018 294

Zubehör:

Gas safety tubing, rubber / Deutsch & Neumann

Safety gas hoses, rubber, without reinforcement / Deutsch & Neumann

Gas lighter / Carl Friedrich Usbeck KG

Piezoelectric gas lighter, Clipper / Deutsch & Neumann

Accessories for Bunsen and Teclu burners / Carl Friedrich Usbeck KG