

Funnel holders, PP

Single or double holder for funnels from 25mm (with adapter supplied) up to 170mm diameter. Without adapter they can be used for powder funnels with stem diameters not exceeding 35mm. Two clamp adapters are provided for use with rods from 8mm to 14mm.

No. of funnels	Min. funnel Ø mm	Max. funnel Ø mm	Funnel	Material	PK	Cat. No.
1	25	170	25 ... 170 mm	PP	1	9.251 561
2	25	170	25 ... 170 mm	PP	1	9.251 560

Alternativen:

Filter funnel stands / BRAND

Funnel stands LaboPlast[®], PP / Bürkle