

Laboratory trays LaboPlast[®], PP

White multi-purpose trays. Injection-moulded. Durable and resist flexing, with strong, rounded rim. Stackable and hygienic. Surface is easy-to-clean. Resistant to many acids and alkaline solutions. Food-safe. Temperature resistant up to 60 °C.

Temperature resistance:

60 °C

Colour:

White

Material:

PP

| Nominal capacity l | Ext. width mm | Ext. depth mm | Ext. height mm | Ext. dim. (WxDxH) | Int. width mm | Depth int. mm | Height int. mm | Int. D (WxD) |
|-----------------------|------------------|------------------|-------------------|----------------------|------------------|------------------|-------------------|-----------------|
| 1.0 | 180 | 225 | 44 | 180 x 225 x 44 mm | 127 | 176 | 39 | 127 x 176 mm |
| 3.0 | 245 | 305 | 69 | 245 x 305 x 69 mm | 175 | 235 | 64 | 175 x 235 mm |
| 5.0 | 305 | 370 | 78 | 305 x 370 x 78 mm | 230 | 295 | 72 | 230 x 295 mm |
| 15.0 | 415 | 515 | 120 | 415 x 515 x 120 mm | 290 | 390 | 117 | 290 x 390 mm |
| 30.0 | 530 | 625 | 135 | 530 x 625 x 135 mm | 410 | 505 | 132 | 410 x 505 mm |
| 60.0 | 645 | 845 | 158 | 645 x 845 x 158 mm | 503 | 710 | 156 | 503 x 710 mm |

Alternativen:

Photographic trays, 18/10 steel / BOCHEM Instrumente

Photographic trays, PVC / Kartell S.p.A.

Laboratory tray LaboPlast[®], set small, PP / Bürkle