

Syringe filter CHROMAFIL[®], Polyamide (PA)

Hydrophilic membrane. The membrane is recommended for filtration of medium polar organic/aqueous liquids.

Sterile:

Connetion inlet:

No
Luer-Lock female

Type	Membrane Ø mm	Pore size µm	Material membrane	Material housing	Autoclavable	Connetion outlet	Colour
CHROMAFIL [®] ; Xtra	13	0.20	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ; Xtra	13	0.45	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ; Xtra	25	0.20	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ; Xtra	25	0.20	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ; Xtra	25	0.45	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ; Xtra	25	0.45	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ;	3	0.20	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ;	3	0.45	PA	PP	Yes	Luer-Slip male	Colourle
CHROMAFIL [®] ;	15	0.20	PA	PP	Yes	Minispike	Yellow/g
CHROMAFIL [®] ;	15	0.45	PA	PP	Yes	Minispike	Green/tr
CHROMAFIL [®] ;	25	0.20	PA	PP	Yes	Luer-Slip male	Yellow/g
CHROMAFIL [®] ;	25	0.20	PA	PP	Yes	Luer-Slip male	Yellow/g
CHROMAFIL [®] ;	25	0.45	PA	PP	Yes	Luer-Slip male	Green/tr
CHROMAFIL [®] ;	25	0.45	PA	PP	Yes	Luer-Slip male	Green/tr