

# LLG-Crimp Neck Micro-Vials ND8

With labelling field:

No

| Nominal capacity<br>ml | Cap size | Ø<br>mm | Height<br>mm | Dimensions (ØxH) | Colour | Bottom shape | Material | PK  | Cat. No.  |
|------------------------|----------|---------|--------------|------------------|--------|--------------|----------|-----|-----------|
| 0.2                    | ND8      | 5.5     | 31.5         | 5.5 x 31.5 mm    | Clear  | Conical      | Glass    | 100 | 6.235 606 |
| 0.3                    | ND8      | 5.5     | 31.5         | 5.5 x 31.5 mm    | Clear  | Round        | Glass    | 100 | 4.001 554 |
| 0.8                    | ND8      | 8.2     | 30.0         | 8.2 x 30 mm      | Clear  | Flat         | Glass    | 100 | 9.003 427 |
| 0.7                    | ND8      | 7.0     | 40.0         | 7 x 40 mm        | Clear  | Conical      | Glass    | 100 | 4.001 515 |
| 0.7                    | ND8      | 7.0     | 40.0         | 7 x 40 mm        | Brown  | Conical      | Glass    | 100 | 6.902 044 |
| 0.7                    | ND8      | 7.0     | 40.0         | 7 x 40 mm        | Clear  | Flat         | Glass    | 100 | 4.008 202 |
| 0.7                    | ND8      | 7.0     | 40.0         | 7 x 40 mm        | Brown  | Flat         | Glass    | 100 | 6.495 247 |
| 1.2                    | ND8      | 8.2     | 40.0         | 8.2 x 40 mm      | Clear  | Flat         | Glass    | 100 | 4.008 206 |
| 1.2                    | ND8      | 8.2     | 40.0         | 8.2 x 40 mm      | Brown  | Flat         | Glass    | 100 | 6.266 864 |

Zubehör:

LLG-Crimper / LLG Labware

LLG-Aluminium Crimp Caps ND8, ready assembled / LLG Labware

LLG-Storage Boxes, PP, with alphanumeric coding / LLG Labware

Alternativen:

LLG-2in1 KITS with Screw Neck Vials ND8 (small opening) / LLG Labware