

Laboratory tubing Isoflex, PVC

Autoclavable, transparent and non-toxic, non ageing and non oxidising. Packed in PE bags. The tubing is marked every metre. Good resistance against strong acids and alkalis, fair resistance against organic solvents. Durometer 77 Shore A.

Colour: Transparent
Material: PVC

Int. Ø mm	Ext. Ø mm	Wall thickness mm	Length m	Shore hardness	PK	Cat. No.
0.8	2.4	0.8	20	77 A	20	9.205 142
2	4	1.0	20	77 A	20	9.205 143
6	9	1.5	20	77 A	20	9.205 146
4	6	1.0	20	77 A	20	9.205 144
5	8	1.5	20	77 A	20	9.205 145
7	11	2.0	20	77 A	20	9.205 147
8	12	2.0	20	77 A	20	9.205 148
10	14	2.0	20	77 A	20	9.205 149
12	17	2.5	20	77 A	20	9.205 150