

behrotest^{®} HTI 1 manual titration station

The HTI 1 manual titration station consists of

-a burette with digital display and

-a magnetic stirrer with a precisely fitting holder for COD reaction vessels. A screen serves as a neutral background and allows the user to determine the colour change at the end of the titration precisely. Therefore a titration could always be carried out under similar visual conditions. This improves the precision and the reproducibility of the results.

The precise positioning of the reaction vessels in the holder on the top of the magnetic stirrer also contributes to this. The angled wing of the screen protects against lateral dazzling light.

Weight: 3.5 kg

Power supply: 115 ... 230 V·AC, 50/60 Hz

Type	Width mm	Depth mm	Height mm	Dimensions (WxDxH)	PK	Cat. No.
HTI	330	200	600	330 x 200 x 600 mm	1	9.920 730