

# Tubing clips, Hoffmann pattern

Brass nickel plated. Bottom hinged.

Int. width mm	Material	PK	Cat. No.
12	Brass-nickel alloy	1	9.180 112
17	Brass-nickel alloy	1	9.180 117
20	Brass-nickel alloy	1	9.180 120
25	Brass-nickel alloy	1	9.180 125
30	Brass-nickel alloy	1	9.180 130
50	Brass-nickel alloy	1	9.180 142

Alternativen:

Tubing clips, Mohr pattern / Carl Friedrich Usbeck KG

Tubing clips, Hoffmann pattern / BOCHEM Instrumente