

Beakers, PP, low form, blue graduations

Highly transparent. Autoclavable to 121 °C for 20 min. and can be used continuously at 100°C. Excellent chemical resistance. Tolerance $\pm 10\%$. Conform to ISO 7056 - 1981 (E).

Nominal capacity ml	Ø mm	Height mm	Standard	Shape	With rim	With spout	With handle	Material	Graduation ml	Type of g
25	34	49	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	1	Printed
50	42	60	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	2	Printed
100	52	72	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	5	Printed
250	71	95	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	10	Printed
500	88	119	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	10	Printed
1000	110	146	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	20	Printed
2000	133	184	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	50	Printed
3000	159	201	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	500	Printed
5000	190	228	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PP	500	Printed

Alternativen:

Beakers, PP, low form, embossed graduations / Kartell S.p.A.