

Flexible Distributors, PP

For connecting elastic tubing by simply slipping it onto hose connectors. Suitable for GL45 threaded bottles. With screw cap made of glass-fibre reinforced PP, ideal for transferring liquids or gases. The angled tubing connectors can be attached to a liquid source, a pressure source or a vacuum system. Avoid crimping of the tubing. The distributor can be loosened by simply turning the screw cap. The complete distribution system can then be fixed onto another bottle.

Colour:

Green

Thread	Bore mm	No. of nozzles	Olives ext. Ø mm	Material	PK	Cat. No.
GL45	6	2	8.8	PP	1	7.634 892
GL45	6	3	8.8	PP	1	9.209 858
GL45	7	4	10.8	PP	1	9.209 859