

Air valves, PP

Red/yellow. With connecting thread UNF $\frac{1}{4}$ " 28G. Filter membrane made of PTFE/PP with a porosity of 1 μm and integrated non-return valve made of silicone. With service life indicator.

Thread	Shelf life	Material	PK	Cat. No.
UNF $\frac{1}{4}$ " 28G	6 months	PP	2	6.313 842
UNF $\frac{1}{4}$ " 28G	6 months	PP	10	6.313 843
UNF $\frac{1}{4}$ " 28G	6 months	PP	50	6.313 844