

# Flanged Tubing, PTFE, PA

Flanged PTFE tubing with black tube end fittings UNF  $\frac{1}{4}$ " 28 G made of PP and washers made of PA. The tubing is ready for use. Temperature resistant from 0 to +100 °C, pressure resistant up to 30 bar, universal chemical resistance.

Material: PTFE  
Operating pressure bar: 30 bar

Int. Ø mm	Ext. Ø mm	Length mm	Max. temp. resistance °C	Temperature resistance	PK	Cat. No.
0.8	1.6	100	100	0 ... 100 °C	1	6.231 140
0.8	1.6	250	100	0 ... 100 °C	1	6.078 162
0.8	1.6	750	100	0 ... 100 °C	1	6.263 025
0.8	1.6	1000	100	0 ... 100 °C	1	6.801 712
1.6	3.2	100	100	0 ... 100 °C	1	6.262 704
1.6	3.2	250	100	0 ... 100 °C	1	6.237 591
1.6	3.2	500	100	0 ... 100 °C	1	6.071 880
1.6	3.2	750	100	0 ... 100 °C	1	6.258 092