

Measuring cylinders, SAN, tall form, class B, moulded graduations

Tall form. Crystal clear. With a raised scale and ring marks at the primary scale points, calibrated IN. Class B tolerances according to DIN 12681/ISO 6706. Hexagonal base with bottom studs provides high stability. Thermal stress up to 60 °C does not cause tolerance limits to be permanently exceeded. Food-safe product according to regulation (EU) No. 10/2011.

Temperature resistance:
Standard:

60 °C
DIN 12681 / ISO 6706

Nominal capacity ml	Division ml	Tolerance ml	Ø mm	Height mm	Max. temp. resistance °C	Shape	Type of graduation	Material	Colour
50	1.0	1.00	28	199	60	Tall	Embossed	SAN	Clear
100	1.0	1.00	34	260	60	Tall	Embossed	SAN	Clear
250	2.0	2.00	47	315	60	Tall	Embossed	SAN	Clear
500	5.0	5.00	61	350	60	Tall	Embossed	SAN	Clear
1000	10.0	10.00	76	415	60	Tall	Embossed	SAN	Clear

Alternativen:

Measuring cylinders, PP, tall form, class B, blue moulded graduations / VITLAB

Measuring cylinders Nalgene™, PMP / Thermo Elect.LED

Measuring cylinders Nalgene™, PMP / Thermo Elect.LED

Graduated cylinders, PP, class B, embossed scale / Kartell S.p.A.