

Syringe filter CHROMAFIL^{®} Xtra, Polyvinylidenfluoride (PVDF)

Hydrophilic membrane, for aqueous solutions, water-soluble oligomers and polymers like proteins binding capacity for proteins 20 µg/25 mm filter.

Connetion inlet: Luer-Lock female

Connetion outlet: Luer-Slip male

Colour: Colourless

Membrane Ø mm	Pore size µm	Material membrane	Material housing	Sterile	Feature	PK	Cat. No.
13	0.20	PVDF	PP	No	Hydrophilic	100	6.287 313
13	0.45	PVDF	PP	No	Hydrophilic	100	6.287 314
25	0.20	PVDF	PP	No	Hydrophilic	100	4.003 413
25	0.20	PVDF	PP	No	Hydrophilic	400	4.003 414
25	0.45	PVDF	PP	No	Hydrophilic	100	9.049 063
25	0.45	PVDF	PP	No	Hydrophilic	400	4.003 415