

Qualitative filter paper, Grade 2294, circles, with central hole

Thick filter card. Fast.

For technical application which range from the filtration of beverages to the purification of electroplating baths.

- Material: Cellulose
- High wet strength
- High flow rate
- Retains medium to coarse particles
- Central hole 33 mm

Filtration speed: 27
Particle retention: 8 ...15 µm

Type	Ø mm	Weight g/m ²	Thickness mm	Max. filtration speed s	Material	Type of analysis	Wet strength	Shape	PK	Cat. No.
2294	180.0	556	1.50	27.5	Cellulose	Qualitative	high	Round filters	100	6.253 836