

# Quantitative filter paper, Grade 589/2, folded filters 1/2

**Medium fast, 4 to 12  $\mu$ m. White ribbon.**

- For the analysis of foodstuffs
- Weight: 85 g/m<sup>2</sup>;
- Thickness: 0.18mm

Ø mm	Weight g/m <sup>2</sup>	Thickness mm	Type of analysis	Shape	PK	Cat. No.
110.0	85	0.18	Quantitative	Folded filters	100	9.060 215
150.0	85	0.18	Quantitative	Folded filters	100	9.060 216

Alternativen:

- Filter paper, qualitative, type MN 615, circles / Macherey-Nagel
- Filter paper, quantitative, round filters, 589/3 / Hahnemühle FineArt