

Support clamps, zinc die casting

Powder-coated. With bosshead. Jaws with cork.
With bosshead:

Yes

Min. clamping width mm	Max. clamping width mm	Clamping range	Material	Type of coating	PK	Cat. No.
12	25	12 ... 25 mm	Die-cast zinc	Powder-coated	1	6.073 337
20	40	20 ... 40 mm	Die-cast zinc	Powder-coated	1	6.053 908
38	60	38 ... 60 mm	Die-cast zinc	Powder-coated	1	6.077 908