

Nuclepore^{®} membrane filter, PC

- Hydrophilic, low protein binding
- Smooth surface for easy identification of particles
- Highly resistant to chemicals and good thermal stability
- Autoclavable up to 121 °C

Applications: Ideal for environmental and water analysis.

Ø mm	Pore size µm	Material	Autoclavable	PK	Cat. No.
25	0.20	PC	Yes	100	6.315 879
25	0.40	PC	Yes	100	6.322 328
25	0.80	PC	Yes	100	6.314 452
25	3.00	PC	Yes	100	6.322 268
25	5.00	PC	Yes	100	6.315 338
25	8.00	PC	Yes	100	6.315 321
25	10.00	PC	Yes	100	6.315 571
25	12.00	PC	Yes	100	4.681 438
47	0.20	PC	Yes	100	6.312 964
47	0.40	PC	Yes	100	6.313 989
47	3.00	PC	Yes	100	6.490 684
50	0.45	PC	Yes	100	6.086 620
50	2	PC	Yes	100	4.684 094
50	5.00	PC	Yes	100	6.312 380
142	1.00	PC	Yes	25	4.681 984