

# Syringe filter CHROMAFIL<sup>®</sup>; Hydrophilized Polytetrafluoroethylene (H-PTFE)

Hydrophilized polytetrafluoroethylene syringe filters labelled for method validation and certification

- Hydrophobic membrane with additional hydrophilic characteristics
- For polar and nonpolar solutions
- Resistant against all kinds of solvents as well as acids and bases

Sterile: No  
Connetion inlet: Luer-Lock female  
Connetion outlet: Luer-Slip male  
Colour: Colourless

| Type                          | Membrane Ø<br>mm | Pore size<br>µm | Material membrane | Material housing | Autoclavable | Feature     | With filter | P  |
|-------------------------------|------------------|-----------------|-------------------|------------------|--------------|-------------|-------------|----|
| CHROMAFIL <sup>®</sup> ; Xtra | 13               | 0.20            | PTFE              | PP               | Yes          | Hydrophilic | Yes         | 10 |
| CHROMAFIL <sup>®</sup> ; Xtra | 13               | 0.45            | PTFE              | PP               | Yes          | Hydrophilic | Yes         | 10 |
| CHROMAFIL <sup>®</sup> ; Xtra | 25               | 0.20            | PTFE              | PP               | Yes          | Hydrophilic | Yes         | 10 |
| CHROMAFIL <sup>®</sup> ; Xtra | 25               | 0.45            | PTFE              | PP               | Yes          | Hydrophilic | Yes         | 10 |
| CHROMAFIL <sup>®</sup> ; Xtra | 25               | 0.45            | PTFE              | PP               | Yes          | Hydrophilic | Yes         | 40 |