

Coloured wash-bottles, narrow-neck, LDPE, set

The set includes one bottle each in red, yellow, blue and green. Screw cap and PP spray tube. Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube.

Colour:

Colour sorted

Nominal capacity ml	Ø mm	Height mm	Dimensions (ØxH)	Thread	Material	Material closure	PK	Cat. No.
500	74	180	74 x 180 mm	GL 25	LDPE	PP	1	7.910 783
1000	92	221	92 x 221 mm	GL 32	LDPE	PP	1	4.671 822

Zubehör:

Wash bottle caps VENT-CAP, PP / VITLAB