

# Magnetic stirrer, operated by water/air pressure

PP-ABS construction. This unique stirrer is operated either by water pressure between 1 kg/cm<sup>2</sup> and 2.5 kg/cm<sup>2</sup>; (15 psi to 40 psi) or by air pressure. No electrical supply is needed, thus eliminating heat, electric shock and explosion hazards. Corrosion-proof, all-plastic construction. Can be used on the bench or clamped to a stand with a 12 mm ( $\frac{1}{2}$ " diameter rod using the integral bosshead.

Dimensions (ØxH): 129 x 61 mm

No. of stirring positions	Material stirring surface	Material housing	Colour housing	Ø mm	Height mm	With timer	PK	Cat. No.
1	PP	ABS	Grey	129	61	No	1	9.303 043