

Beakers, PMP (TPX[®]), low form

Crystal clear. Conforms to ISO 7056 - 1981 (E). Autoclavable. Will withstand 170°C for short periods. Permanent blue printed graduations. Excellent chemical resistance. Tolerance \pm 10%.

Nominal capacity ml	Ø mm	Height mm	Standard	Shape	With rim	With spout	With handle	Material	Graduation ml	Type of g
25	34	49	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	1	Printed
50	42	60	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	2	Printed
100	52	72	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	5	Printed
250	71	95	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	10	Printed
500	88	119	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	10	Printed
1000	110	146	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	20	Printed
2000	133	184	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	50	Printed
3000	159	201	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	500	Printed
5000	190	228	ISO 7056 - 1981 (E)	Low	Yes	Yes	No	PMP	500	Printed