

# Gravimetric dilutor DiluFlow<sup>&reg;</sup>, with double pump

All stainless steel, ultra-low profile for easy use under laminar flow.

- Quick and accurate
- Compact and user-friendly
- RS232 output for LIMS and printer, USB

**Scope of supply:** 1 BagOpen<sup>&reg;</sup> 400 DiluFlow<sup>&reg;</sup>,  
1 fixed arm, 1 power cord, 1 &oslash; 6.4 mm tubings set per pump for GL45 bottle, 1 drip  
tray, 1 user's manual, guarantee

Material housing:	Stainless steel
Dimensions (WxDxH):	370 x 437 x 394 mm
Weight:	12.5 kg
Power supply:	100 ... 240 V, 50/60 Hz
IP code:	IP 22
Warranty:	1 year
Min. voltage:	100 V
Max. voltage:	240 V
Power supply frequency:	50/60 Hz

Type	Dilution factor	Resolution g	Min. weighing range g	Max. weighing range g	Width mm	Depth mm	Height mm	Type of auxiliary ener
DiluFlow&reg;	1:2 ... 1:99	0.01	0.1	3000	370	437	394	Mains connect